DECIFERA

FOAM SILICONE BORDER

EXTRA THIN

EN

User manual for silicone foam dressing with border, extra thin

Description

Silicone foam dressing is an absorbent, atraumatic dressing made from polyurethane foam. The outer surface of the foam is bonded to a vapour-permeable polyurethane membrane, which acts as a barrier to liquid and microorganisms. The wound contact surface of Silicone foam dressing is coated with a layer of soft silicone that does not stick to the surface of a wound or cause trauma to delicate new tissue upon removal. This soft silicone layer is also slightly tacky, which facilitates application and retention of the dressing to intact skin, but does not cause epidermal stripping or pain on removal. The centre of the dressing consists of polyurethane foam that can lock exudate into the product, preventing the wound bed from drying out.

Indications

The product is suitable for many types of exuding wounds, including leg and pressure ulcers, superficial and partial thickness burns, donor sites, postoperative wounds and skin abrasions. The dressing absorbs exudate and maintains a moist wound-healing environment whilst minimising the risk of maceration.

Contra-indications

Discontinue use if the patient is allergic to any components of the dressing.

Application method

1. Make sure the wound area is clean and dry before using the dressing.
2. A dressing should be selected that overlaps the wound margin by at least two centimetres.
3. Remove the release film before use and adhesive the dressing on the wound area.

Frequency of dressing changes

The interval between changes will normally be determined by the degree of exudate produced. The dressing may be left up to 7 days.

Warnings

Sloughy wounds dressed with silicone foam dressing may initially appear to increase in size due to autolytic debridement promoted by the moist conditions produced beneath the dressing. This is normal and to be expected. Not for use with oxidizing solutions, such as Hypochlorite solution or Hydrogen peroxide solution (H₂O₂).

Storage

Store at a normal room temperature, 5 - 35 °C, in a well-ventilated clean room, with a relative humidity of < 80% and no corrosive gas.

Other

  Silicone foam dressings have been sterilized using EO.

-  Do not reuse or resterilize product.
-  Do not use if package has been compromised.
-  Protect against direct sunlight.
-  Keep dry.
-  Store at a normal room temperature 5 - 35 °C.

DECIFERA® is a registered trademark of Van Heek Medical. All rights reserved.

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